

## COLD HYDROGEL FOR FEVER



### PRODUCT DESCRIPTION

**Cold Hydrogel for fever** is a hydrogel patch that provides cooling relief effect to children when they have slight fever.

Patches will stay cool for up to 8 hours and will remain in place on the skin without leaving a sticky residue. It is comfortable and easy to apply.

### INDICATIONS

The product is designed for use as an adjunctive in the case of fever, providing relief through a pleasant cooling effect, thereby being helpful in common antipyretic therapies.

The local reduction in the temperature of the forehead improves conditions for rest.

### MECHANISM OF ACTION

The contact of the hydrogel layer with the skin allows the device itself to exert cryogenic action. This is due to the slow evaporation of the water contained in the gel that, because of heat absorption, determines the reduction in temperature of the forehead.

### INSTRUCTIONS FOR USE

1) Clean the affected area prior to application to not hinder adhesion.

2) Do not apply to wounds or damaged skin.

### PRESENTATION

1 printed foil sachet containing 2 patches.

Patch size: 40 mm x 60 mm.

Sachet size: 80 mm x 130 mm.

### PRODUCT CATEGORY

Medical Device I.